



Product Datasheet

SGK1 (Phospho-S78) antibody (orb393161)

Description	Rabbit polyclonal antibody to SGK1 (Phospho-S78)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S78 of human SGK1 protein. The exact sequence is proprietary.
Target	SGK1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q06226 , Q9WVC6 , O00141
Entrez	20393 , 6446
Dilution Range	WB: 1:500:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of SGK1 (Phospho-S...



Immunohistochemical analysis of SGK1 (Ph...