



Product Datasheet

HNF1A antibody (orb39410)

| | |
|----------------------------|--|
| Description | Purified rabbit polyclonal antibody. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IF, WB |
| Immunogen | Synthetic Peptide |
| Target | HNF1A |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 67356 |
| Uniprot ID | P20823 |
| NCBI | NP_000536.5 |
| Dilution Range | WB: 1:1000, IF/ICC: 1:10-50 |
| Expiration Date | 12 months from date of receipt. |



Confocal immunofluorescent analysis of H...