


## Product Datasheet

### 6\*His, His-Tag antibody (orb394654)

<b>Description</b>	Mouse monoclonal antibody to 6*His, His-Tag.
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IP, WB
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>Dilution Range</b>	WB: 1:1000-1:10000, IP: 1:1000-1:10000, IF/ICC: 1:20-1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of 6\*His-tagged fu...