


Product Datasheet

6*His, His-Tag antibody (orb394655)

Description	Mouse monoclonal antibody to 6*His, His-Tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 6*His-tagged fu...