


Product Datasheet

6*his,His tag antibody (orb394656)

Description	Rabbit polyclonal antibody to 6*his,His tag.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of 6*His-tagged fu...