



Product Datasheet

Alpha galactosidase A antibody (orb394691)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Alpha galactosidase A. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, IP, WB |
| Immunogen | galactosidase, alpha |
| Target | GLA |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 49 kDa |
| Uniprot ID | P06280 |
| Dilution Range | WB:1:500-1:5000, IHC:1:20-1:200, IP:1:500-1:5000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HeLa cells using...



Immunohistochemical staining of paraffin...