



Product Datasheet

Coproporphyrinogen oxidase antibody (orb394905)

Description	Rabbit polyclonal antibody to Coproporphyrinogen oxidase.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	coproporphyrinogen oxidase
Target	CPOX
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	36 kDa
Uniprot ID	P36551
Dilution Range	WB: 1:1000-1:10000 IP: 1:200-1:2000 IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human liver tis...



Immunohistochemical staining of paraffin...