



Product Datasheet

IL18RAP antibody (orb39492)


Description	Rabbit polyclonal antibody to IL18RAP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IL18RAP
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Purity	This antibody is purified through a protein A column, followed by peptide affinity purification.
MW	68310 Da
Uniprot ID	O95256
NCBI	NP_003844.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of k562 cells
u...



Immunohistochemical staining of paraffin...



Western blot analysis of in K562 cell li...