

## Product Datasheet

### GLS antibody (orb39504)

**Description** Peptide affinity purified rabbit polyclonal antibody.

**Species/Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** FC, IF, IHC-P, WB

**Immunogen** Synthetic Peptide

**Target** GLS

**Form/Appearance** Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

**Note** For research use only

**Isotype** Rabbit IgG

**Clonality** Polyclonal



**MW** 73461

**Uniprot ID** [O94925](#)


**NCBI** [NP\\_055720.3](#), [NP\\_001243239.1](#)

**Dilution Range** WB: 1:1000, IHC-P: 1:10-50, IF/ICC: 1:10-50, FACS: 1:10-50

**Expiration Date** 12 months from date of receipt.



Immunohistochemical analysis of formali...



Flow cytometric analysis of GLS antibody...