


Product Datasheet

GFP tag antibody (orb395079)

Description	Rabbit polyclonal antibody to GFP tag.
Species/Host	Rabbit
Reactivity	Human, Jellyfish
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	GFP tag
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Dilution Range	WB: 1:1000-1:4000; IP: 1:2000-1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of eGFP-tagged fus...