

Product Datasheet

Giantin antibody (orb395080)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Giantin. |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF |
| Immunogen | golgin B1, golgi integral membrane protein |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 376 kDa |
| Uniprot ID | Q14789 |
| Dilution Range | IF: 1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |