



Product Datasheet

GRP94 antibody (orb395101)

Description	Mouse monoclonal antibody to GRP94.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	heat shock protein 90kDa beta(Grp94), member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	94 kDa
Uniprot ID	P14625
Dilution Range	WB: 1:500-1:3000; IHC: 1:50-1:500; IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical staining of paraffin...