

Product Datasheet

ICP27 antibody (orb395163)

Description	Rabbit polyclonal antibody to ICP27.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	ICP27
Target	ICP27
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	55 kDa
Uniprot ID	P10238
Expiration Date	12 months from date of receipt.