


Product Datasheet

Integrin alpha-3 antibody (orb395202)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to Integrin alpha-3. |
| Species/Host | Mouse |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | integrin, alpha 3(antigen CD49C, alpha 3 subunit of VLA-3 receptor) |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 120 kDa |
| Uniprot ID | P26006 |
| Dilution Range | WB: 1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of PC-3 cells usin...