



Product Datasheet

MGEA5 antibody (orb395307)

Description	Mouse monoclonal antibody to MGEA5.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	meningioma expressed antigen 5(hyaluronidase)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	130 kDa
Uniprot ID	O60502
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Human brain lys...



Immunohistochemical staining of paraffin...