



## Product Datasheet

### Myosin Light Chain 2 antibody (orb395332)

<b>Description</b>	Mouse monoclonal antibody to Myosin Light Chain 2.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, WB
<b>Immunogen</b>	myosin, light chain 2, regulatory, cardiac, slow
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>MW</b>	19 kDa
<b>Uniprot ID</b>	<b>P10916</b>
<b>Dilution Range</b>	WB: 1:1000-1:5000; IHC: 1:100-1:1000; IF: 1:50-1:500
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical staining of paraffin...