


Product Datasheet

ORM2-Specific antibody (orb395389)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to ORM2-Specific. |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, WB |
| Immunogen | orosomuroid 2 |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 45 kDa |
| Uniprot ID | P19652 |
| Dilution Range | WB: 1:500-1:5000; IF: 1:50-1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of human plasma ti...