



Product Datasheet

Akt (Phospho-Tyr315) antibody (orb395432)

Description	Mouse monoclonal antibody to Akt (Phospho-Tyr315).
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	v-akt murine thymoma viral oncogene homolog 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	60KD
Uniprot ID	P31749
Dilution Range	WB: 1:1000-1:10000; IHC: 1:100-1:400
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...