

Product Datasheet

RCC1 (phospho-Ser12) antibody (orb395435)

Description Rabbit polyclonal antibody to RCC1 (phospho-Ser12).

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen regulator of chromosome condensation 1

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 48 kDa


Uniprot ID **P18754**

Dilution Range WB: 1:500-1:1000; IHC: 1:20-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of HeLa cells using RCC1 (phospho-Ser12) antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using RCC1 (phospho-Ser12) antibody.