



Product Datasheet

PLK1 (phospho-S326) antibody (orb395455)

Description	Rabbit polyclonal antibody to PLK1 (phospho-S326).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	polo-like kinase 1(Drosophila)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	P53350
Dilution Range	WB: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 cells using PLK1 (phospho-S326) antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using PLK1 (phospho-S326) antibody.