



Product Datasheet

POLR2B-Specific antibody (orb395459)

Description	Rabbit polyclonal antibody to POLR2B-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	polymerase(RNA) II(DNA directed) polypeptide B, 140kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P30876
Dilution Range	WB: 1:500-1:5000; IP: 1:200-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells using POLR2B-Specific antibody.



Immunohistochemical staining of paraffin-embedded tissue using POLR2B-Specific antibody.