



Product Datasheet

RABEP2 antibody (orb395501)

Description	Mouse monoclonal antibody to RABEP2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	rabaptin, RAB GTPase binding effector protein 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	62 kDa
Uniprot ID	Q9H5N1
Dilution Range	WB: 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 cells using RABEP2 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using RABEP2 antibody.