


Product Datasheet

SEC61alpha antibody (orb395545)

Description	Rabbit polyclonal antibody to SEC61alpha.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Sec61 alpha 1 subunit(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	52 kDa
Uniprot ID	P61619
Dilution Range	WB: 1:200-1:2000 IF/ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...