

Product Datasheet

SEC61alpha antibody (orb395546)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to SEC61alpha. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, WB |
| Immunogen | Sec61 alpha 2 subunit(S. cerevisiae) |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 52 kDa |
| Uniprot ID | P61619 |
| Dilution Range | WB: 1:500-1:2000; IF: 1:20-1:200 |
| Expiration Date | 12 months from date of receipt. |