


Product Datasheet

STXBP5 antibody (orb395627)

Description	Rabbit polyclonal antibody to STXBP5.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	syntaxin binding protein 5(tomosyn)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	130 kDa
Uniprot ID	Q5T5C0
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...