



Product Datasheet

DNAJA2 antibody (orb39565)


| | |
|----------------------------|---|
| Description | DNAJA2 antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | Recombinant Human Dnaj homolog subfamily A member 2 protein (233-412AA) |
| Target | DNAJA2 |
| Preservatives | PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | O60884 |
| Dilution Range | WB:1:1000-5000, IHC-P:1:20-200 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of A375(lane 1), P...



Immunohistochemical staining of human CO...



Immunohistochemical staining of human to...