



## Product Datasheet

### TFIP11 antibody (orb395671)

<b>Description</b>	Mouse monoclonal antibody to TFIP11.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	tuftelin interacting protein 11
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>MW</b>	97 kDa
<b>Uniprot ID</b>	<b>Q9UBB9</b>
<b>Dilution Range</b>	WB: 1:200-1:1000; IHC: 1:50-1:500
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa cells usin...



Immunohistochemical staining of paraffin...