


Product Datasheet

TTC7B antibody (orb395726)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TTC7B. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | tetratricopeptide repeat domain 7B |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 96 kDa |
| Uniprot ID | Q86TV6 |
| Dilution Range | WB: 1:200-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of mouse liver tis...