



## Product Datasheet

### TXNL6 antibody (orb395734)

<b>Description</b>	Rabbit polyclonal antibody to TXNL6.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, IP, WB
<b>Immunogen</b>	nucleoredoxin-like 1
<b>Target</b>	NXNL1
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>MW</b>	17 kDa
<b>Uniprot ID</b>	<a href="#">Q96CM4</a>
<b>Dilution Range</b>	WB:1:500-1:5000, IHC:1:20-1:200, IP:1:200-1:2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Immunohistochemical staining of paraffin...