

Product Datasheet

UL27 antibody (orb395745)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to UL27. |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Envelope glycoprotein B |
| Target | gB |
| Form/Appearance | liquid |
| Storage | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.) |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | ≥95% as determined by SDS-PAGE |
| MW | 80 kDa |
| Uniprot ID | P10211 |
| Expiration Date | 12 months from date of receipt. |