

Product Datasheet

Vpr antibody (orb395761)

Description	Rabbit polyclonal antibody to Vpr.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	vpr
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	33 kDa
Uniprot ID	P05928
Expiration Date	12 months from date of receipt.