



Product Datasheet

WNT5A-B antibody (orb395771)

Description	Rabbit polyclonal antibody to WNT5A-B.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	wingless-type MMTV integration site family, member 5A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	42 kDa
Uniprot ID	P41221
Dilution Range	WB:1:200-1:2000, IHC:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells usi...



Immunohistochemical staining of paraffin...