



Product Datasheet

CD20 antibody (orb396140)


Description	Mouse monoclonal antibody to CD20.
Species/Host	Human
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human lymphoblastoid cell line SB.
Target	CD20
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P11836
Expiration Date	12 months from date of receipt.



Flow-cytometry analysis of Rhesus monkey...



Flow-cytometry analysis of Cynomolgus mo...



Flow-cytometry analysis of Human lymphoc...