

## Product Datasheet

### CD133E1 antibody (orb397155)

<b>Description</b>	CD133E1 antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Immunogen</b>	Protein
<b>Target</b>	CD133E1
<b>Form/Appearance</b>	10 mM PBS, 0.02% sodium azide
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purified by Protein A
<b>Expiration Date</b>	12 months from date of receipt.