

Product Datasheet

phospho-MLKL (S358) Mouse Monoclonal Antibody(6F8) (orb397214)

Description	phospho-MLKL (S358) Mouse Monoclonal Antibody(6F8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q8NB16
Expiration Date	12 months from date of receipt.