




Product Datasheet

ANGPTL4 antibody (orb39902)

Description	Rabbit polyclonal antibody to ANGPTL4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ANGPTL4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	45214
Uniprot ID	Q9BY76
NCBI	NP_647475.1 , NP_001034756.1
Dilution Range	WB: 1:2000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formali...



Flow cytometric analysis of ANGPTL4 anti...