



Product Datasheet

AAT Polyclonal Antibody (orb399394)


Description	AAT Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P01009
Expiration Date	12 months from date of receipt.



Western Blot analysis of NIH-3T3 cells u...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...