

Product Datasheet

CD11c antibody (FITC) (orb399682)

| | |
|----------------------------|--|
| Description | CD11c antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | FITC |
| Tested Applications | ICC, IF |
| Immunogen | KLH conjugated synthetic peptide derived from human CD11c. |
| Target | CD11c |
| Form/Appearance | Liquid: PBS, 2 mM Sodium Azide, 30% glycerol |
| Concentration | - 100 µg (in 200 µl): 0.5 mg/ml |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Polyclonal antibodies are purified by peptide affinity chromatography |
| MW | 128 kDa |
| Uniprot ID | P20702 , Q9QXH4 |
| Entrez | 3687 |
| Dilution Range | IF/ICC: 1:100-300 |
| Expiration Date | 12 months from date of receipt. |