
Product Datasheet

CD169_Sialoadhesin antibody (FITC) (orb399691)

Analysis Reagents

CD169_Sialoadhesin antibody (FITC)

Description	CD169_Sialoadhesin antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CD169_Sialoadhesin.
Target	CD169_Sialoadhesin (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	183 kDa
Uniprot ID	Q62230 , Q9BZZ2
NCBI	NP_075556.1 , NM_023068.3
Entrez	6614
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.