
Product Datasheet

CD200 antibody (FITC) (orb399694)

Description

CD200 antibody (FITC)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthetic peptide derived from human CD200.

Target

CD200 (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

31 kDa

Uniprot ID[P41217](#)**NCBI**[NM_001004196.2](#), [NP_001004196.21](#), [NM_001004196.21](#), [NP_001004196.2](#)**Entrez**[4345](#)**Dilution Range**

IF/ICC: 1:100-300

Expiration Date

12 months from date of receipt.