
Product Datasheet

CD69 antibody (FITC) (orb399709)

Analysis Reagents.

CD69 antibody (FITC)

Description	CD69 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human CD69.
Target	CD69 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	23 kDa
Uniprot ID	Q07108 , P37217
NCBI	NP_001772.1 , NM_001781.21 , NM_001781.2 , NP_001772.11
Entrez	969
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.