

Product Datasheet

IL26 antibody (FITC) (orb399812)

Description	IL26 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human IL26.
Target	IL26 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	20 kDa
Uniprot ID	Q9NPH9
Entrez	55801
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.