

## Product Datasheet

### MASP2 antibody (FITC) (orb399831)

<b>Description</b>	MASP2 antibody (FITC)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human MASP2.
<b>Target</b>	MASP2 (FITC)
<b>Form/Appearance</b>	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Polyclonal antibodies are purified by peptide affinity chromatography
<b>MW</b>	76 kDa
<b>Uniprot ID</b>	<a href="#">O00187</a> , <a href="#">Q9JJS8</a> , <a href="#">Q91WPO</a>
<b>NCBI</b>	<a href="#">NM_006610.3</a> , <a href="#">NP_006601.2</a>
<b>Entrez</b>	<a href="#">10747</a>
<b>Dilution Range</b>	IF/ICC: 1:100-300
<b>Expiration Date</b>	12 months from date of receipt.