

Product Datasheet

Reg3 gamma antibody (FITC) (orb399904)

Description

Reg3 gamma antibody (FITC)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Reg3 gamma.
Target	Reg3 gamma (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	19 kDa
Uniprot ID	O09049 , Q6UW15
NCBI	NM_001008387.2 , NM_001008387.21 , NP_001008388.1 , NP_001008388.11
Entrez	130120
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.