


Product Datasheet

hist2h2l antibody (orb400300)

Description	hist2h2l antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish Histone H2B 3 protein (2-126AA)
Target	hist2h2l
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PC60
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Zebrafish tissue...