

Product Datasheet

hist2h2l antibody (HRP) (orb400301)

Description	hist2h2l antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Histone H2B 3 protein (2-126AA)
Target	hist2h2l
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PC60
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.