



## Product Datasheet

### COPS3 antibody (orb40031)


<b>Description</b>	COPS3 antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, IP, WB
<b>Immunogen</b>	Recombinant Human COP9 signalosome complex subunit 3 protein (194-423AA)
<b>Target</b>	COPS3
<b>Preservatives</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q9UNS2</a>
<b>Dilution Range</b>	WB:1:1000-5000, IHC-P:1:20-200, IP:1:200-2000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), 2...



Immunohistochemical staining of human te...



Immunoprecipitation analysis of HeLa who...