

## Product Datasheet

### TRIP10 antibody (FITC) (orb400658)

|                        |   |
|------------------------|---|
| <b>Description</b>     | TRIP10 antibody (FITC)  |
| <b>Species/Host</b>    | Rabbit  |
| <b>Reactivity</b>      | Human   |
| <b>Conjugation</b>     | FITC  |
| <b>Immunogen</b>       | Recombinant Human Cdc42-interacting protein 4 protein (276-486AA) |
| <b>Target</b>          | TRIP10  |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol    |
| <b>Form/Appearance</b> | Liquid  |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.     |
| <b>Note</b>            | For research use only   |
| <b>Isotype</b>         | IgG   |
| <b>Clonality</b>       | Polyclonal  |
| <b>Uniprot ID</b>      | <a href="#">Q15642</a>  |
| <b>Dilution Range</b>  | ELISA   |
| <b>Expiration Date</b> | 12 months from date of receipt.                                   |