



Product Datasheet

NDUFS4 antibody (orb400888)


Description	NDUFS4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] iron-sulfur protein 4, mitochondrial protein (43-161AA)
Target	NDUFS4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43181
Dilution Range	WB:1:500-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of U87 whole cell ...



Immunohistochemical staining of human CO...



Immunohistochemical staining of Human ki...