


Product Datasheet

Pcsk1 antibody (orb400896)

Description	Pcsk1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rat Neuroendocrine convertase 2 protein (109-637AA)
Target	Pcsk2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28841
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...